

BITUMEN MC 3000

Characteristic	Value	Test Method
Kinematic Viscosity at 60 °C (cSt)	3000-6000	ASTM D-2170
Flash Point (top open cup) (°C)	66 min	ASTM D-3143
Distillate, %Volume of Total Distillate to at 360°C:		
To 225°C	---	ASTM D-402
To 260°C	0-15	
To 316°C	15-75	
Residue from Distillation to 360°C (%Vol by difference)	80 min	ASTM D-402
Test on Residue from Distillation		
Viscosity at 60°C	30-120	ASTM D-2170
Penetration at 25°C (dmm)	120-250	ASTM D-5
Ductility at 25°C (cm)	100 min	ASTM D-113
Solubility in Trichloroethylene (%wt)	99 min	ASTM D-2042
Water Content (%vol)	0.2 max	ASTM D-95