

BITUMEN MC 70

Characteristic	Value	Test Method
Kinematic Viscosity@ 60°C (cSt)	70-140	ASTM D-2170
Flash Point (Top open cup) (°C)	38 min	ASTM D-3143
Distillate, %volume of total distillate to @ 360°C: to 225°C to 260°C to 316°C	20 max 20-60 65-90	ASTM D-402
Residue from Distillation to 360°C (%Vol by difference)	55 min	ASTM D-402
Test on Residue from Distillation		
Viscosity@ 60°C (cSt)	30-120	ASTM D-2170
Penetration @ 25°C (dmm)	120-250	ASTM D-5
Ductility @ 25°C (cm)	100 min	ASTM D-113
Solubility in Trichloroethylene (%wt)	99 min	ASTM D-2042
Water Content (%vol)	0.2 max	ASTM D-95